# SMART PRESCRIPTION APPLICATION

## Software Design

By

Mr. Natthakan Kaeokanpai 552115020

Mr. Phithiwat Sitthitun 552115051

Department of Software Engineering
College of Arts, Media and Technology
Chiang Mai University

Project Advisor

\_\_\_\_\_

Dr.Prompong Sugunnasil

## **Document History**

Document Name	Version	Status	Date	Viewable	Editable	Responsible
Documents						
Smart prescription application _Software design	Software Design _V0.1	Reviewed	June 10th 2016	NAT,PHI, PROM	NAT,PHI	NAT,PHI
Smart prescription application _Software design	Software Design _V1.0	Released	July 10th 2016	NAT,PHI, PROM	NAT,PHI	NAT,PHI
Smart prescription application _Software design	Software Design _V.1.1	Released	July 20th 2016	NAT,PHI, PROM	NAT,PHI	NAT,PHI
Smart prescription application _Software design	Software Design _V.2.0	Released	July 30th 2016	NAT,PHI, PROM	NAT,PHI	NAT,PHI

#### Acronym

NAT = Natthakan Kaeokanpai

PHI = Phithiwat Sitthitun

PROM = Prompong Sugunnasil

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	2
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

## Content

Chapter One   Introduction	5
1.1 Objective	5
1.2 Project Scope	5
1.3Acronyms and Definitions	5
1.3.1 Acronyms	5
1.3.2Definitions	6
Chapter Two   System Architecture	7
Chapter Three   Detailed Design	8
3.1 Class Diagram of smart prescription application	8
3.2 Class Description	10
(Mobile application)	27
Chapter Four   Database Design	45
4.1 ER Diagram	45
4.2 Table Description	46
Chapter Five   Sequence Diagram	49
SD-01: Doctors can view a doctor home page on the web application	49
SD-02: Doctors can create a patient's profile on the web application	50
SD-03: Doctors can update a patient's profile on the web application	51
SD-04: Doctors can delete a patient's profile on the web application.	52
SD-05: Doctors can search a patient's profile on the web application.	53
SD-06: Doctors can view a patient's profile on the web application	54
SD-07: Doctors can view a list of patient's profile on the web application	55
SD-08: Administrators can view an admin home page on the web application.	56
SD-09: Administrators can create a user's profile on the web application	57
SD-10: Administrators can update a user's profile on the web application	58
SD-11: Administrators can delete a user's profile on the web application	59
SD-12: Administrators can search a user's profile on the web application.	60

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	3
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

SD-13: Administrators can view a user's profile on the web application	61
SD-14: Administrators can view a list of user's profile on the web application	62
SD-15: Doctors, Administrators, and FDAs can login to the web application	63
SD-16: Doctors, Administrators, and FDAs can logout from the web application	64
SD-17: Patients can login to the mobile application.	65
SD-18: Patients can logout from the mobile application.	66
SD-19: Patients can view the QR code on the mobile application.	67
SD-20: Pharmacists can scan the QR code on the mobile application	68
SD-21: Pharmacists can add the time of dispensation to the patient's profile on the mobile application.	69
SD-22: Doctors can create an allergy report on the web application	70
SD-23: Doctors can update an allergy report on the web application	71
SD-24: Doctors can delete an allergy report on the web application	72
SD-25: Doctors can search an allergy report on the web application.	73
SD-26: Doctors can view an allergy report on the web application.	74
SD-27: Doctors can view a list of allergy report on the web application.	75
SD-28: FDAs can view a FDA home page on the web application.	76
SD-29: FDAs can view a list of allergy report on the web application.	77
SD-30: FDAs can view a list of allergy report on the web application.	78
Chapter Six   User Interface Design	79

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	4
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

## **Chapter One | Introduction**

## 1.1 Objective

The objective of the Software Design Document for "Smart prescription application" both Android mobile application and web application are to provide a description and explain of the design system including activity following requirements specification and before programming. This document is based on the project proposal, project plan and software requirement specification and help developers to understand how to work for this application.

### 1.2 Project Scope

Smart prescription application consists of a web application and a mobile application. Smart prescription application requires a web application for creating a patient's profile and allergy report by a doctor and the Smart prescription application require a mobile application to generate and scan the QR code on a mobile phone by a pharmacist.

The features of the web application consist of:

Feature#1: Prescription system

Feature#2: Account management system

Feature#3: Authentication system Feature#5: Report allergy system

Feature#6: Allergy drug summary report system The features of the mobile application consist of: Feature#4: Prescription verification system

#### 1.3Acronyms and Definitions

#### 1.3.1 Acronyms

SDD Software Design Document

CD Class Diagram

URS User Requirement Specification

UI User Interface

SD Sequence Diagram

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	5
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

#### 1.3.2Definitions

#### Class diagram

A type of static structure diagram that describes the structure of a system by showing the system's classes, their attributes, operations (or methods), and the relationships among objects.[1]

#### **Feature**

Activity after input parameters and output the processing to user interface. It explain about functionality of an application used for requirements analysis, coding, design, and testing

#### Sequence diagram

A sequence diagram is an interaction diagram that shows how processes operate with one another and in what order. A sequence diagram shows object interactions arranged in time sequence.

#### **UML**

A standard notation for the modeling of real-world objects as a first step in developing an object-oriented design methodology [2]

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	6
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

## **Chapter Two | System Architecture**




Figure 1: System Architecture (Overview)

The architecture of "Smart Prescription Application", show in the figure 1, consists of two parts. The first part is the web application. In this system, there are three types of user doctor, FDA and admin. The users can use web application via a web browser by access to the web application by login. Doctors can create the patient's profiles and allergy report to MySQL database. The FDAs can view the allergy report on the web application. An admin can create and manage account of doctor and FDA to MySQL database.


The second part is the mobile application Smart Prescription Application stand for android operating system. There have two types of users; patients and pharmacists. Patients can view the QR code on the mobile application. Pharmacist can verify QR code and get the information of patient's profile by scan QR code function from database on the mobile application.

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	7
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016


## **Chapter Three | Detailed Design**

## 3.1 Class Diagram of smart prescription application

(Web application)



### (Mobile application)



<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	9
<b>Document Type</b>	Software Design	Release Date	July 30th, 2016	Print Date	July 30 <sup>th</sup> , 2016

## 3.2 Class Description

(Web application)

## 3.2.1 Class-01: LoginForm

LoginForm
-\$usemame: string
-\$password: string

#### Description:

The LoginForm model class is using for store the username and password to login.

#### Attributes:

No.	Name	Description	Туре
1	\$username	This attribute is username of users	String
2	\$password	This attribute is password of users	String

#### 3.2.2 Class-02: SiteController

SiteController
+login(String username, String password) : Boolean +logout() : Boolean
Powered By Visual Paradigm Community Edition

#### Description:

This class controls the login and logout function method.

No.	Name	Description	Parameter Lists	Return Type
1	login()	<ul> <li>This method is logs in the user using the given username and password in the model.</li> <li>This method returns true when users input username and password are correct.</li> <li>This method returns false when users input username or password is incorrect.</li> </ul>	String username String password	Boolean
3	logout()	- This is the action to logs out the current user and redirect to homepage.	-	Boolean

#### 3.2.3 Class-03: User

User	
-\$id: int	
-\$name : string	
-\$surname : string	
-\$usemame : string	
-Spassword : string	
-\$password2 : string	
-Sroleld : int	

## Description:

This is the User model class for table "tbl\_user". It used to store the user's information. This model class manages and invokes the user's information.

#### Attributes:

No.	Name	Description	Туре
1	\$id	This attribute is id of user.	Int
2	\$name	This attribute is name of user.	String
3	\$surname	This attribute is surname of user.	String
5	\$username	This attribute is username of user.	String
6	\$password	This attribute is password of user.	String
7	\$password2	This attribute is confirm password of user.	String
8	\$roleId	This attribute is role of user, which include 1, 2, and 3.	Int

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	11
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

## 3.2.4 Class-04: UserController

	UserController
+getUs	serList(): ArrayList
+viewL	Jser(Int userId) : Mixed
	eUser(User user): Boolean
+updat	eUser(User user) : Boolean
+delete	eUser(Int userId) : Boolean
+searc	hUser(String users) : ArrayList
	odel(Int id): Model

## Description:

This class is UserController using for controls and manages (CRUD) the user's profile and information. This class controls the user by admin.

No.	Name	Description	Parameter Lists	Return Type
1	getUserList()	<ul><li>This method is used to providing the list all users.</li><li>This method is used to display over all list of user that exist on the database.</li></ul>	-	ArrayList <user></user>
2	viewUser(Int userId)	<ul> <li>userId is id of user from user table.</li> <li>This method used for view the single user's information which includes id, position, name, surname, username, and password to display on the 'view' page by using id of user.</li> <li>The mechanism of this method involves other method: findModel(int id).</li> <li>This method returns user's information when user id that exist on the user table via user model.</li> <li>This method returns false if the user id that does not exist on the user table.</li> </ul>	Int userId	Mixed

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	12
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016


3	createUser(User user)	- This method is used to creating the new data of user into database, the data that input to create will be validated by AJAX validation.  - The mechanism of this method involves other method: validate() and save() from framework.  - This method returns true if admin input overall information is correct.  -This method returns false if admin does not input some information or inputting is incorrect.	User user	Boolean
4	updateUser(User user)	<ul> <li>This method is used to updating the data of user from the database by using user id, the data that input will be validated by AJAX validation.</li> <li>The mechanism of this method involves other method: fineModel(int id), validate(), and save() from framework.</li> <li>This method returns true if admin input the new overall information is correct.</li> <li>This method returns false if admin does not input some information or inputting is incorrect.</li> </ul>	User user	Boolean

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	13
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

5	deleteUser(Int	- userId is id of user from user table.	Int userId	Mixed
	userId)	- This method is used to deleting the user from the database by using delete() method from framework via user id.		
		- The mechanism of this method involves other method: findModel(int id) and delete() from framework.		
		- This method returns true when delete user id that exist on the user table.		
		- This method will return false when delete user id that does not exist on the user table.		
6	searchUser(String users)	- This method is search query applied.	String users	ArrayList <user></user>
		- This method to get ActiveDataProvider instance which will filter models according to data in model fields.		
		- This method returns data items currently available when input attributes that exist in the database.		
7	findModel(Int id)	- id is user id from user table.	Int id	Model
		- This method is used to finding the model based on its primary key value.		
		- If the model is not found, a 404 HTTP exception will be thrown.		
		- This method returns model object by id.		

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	14
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

#### 3.2.5 Class-05: Role




## Description:

This is the Role model class for table "tbl\_role". This model class used to store the role information.

#### Attributes:

No.	Name	Description	Туре
1	\$id	Id of role	Int
2	\$name	The name of role	String

#### 3.2.6 Class-06: AdminController




#### Description:

The AdminController class is the controller to separates the particular user on role type of admin and lead to the specific index page about the User model.

No.	Name	Description	Parameter Lists	Return Type
1	findUserModel(Int id)	<ul> <li>id is user id from user table.</li> <li>This method is used to finding the model based on its primary key value of User model.</li> <li>If the model is not found, a 404 HTTP exception will be thrown.</li> <li>This method returns model object by id.</li> </ul>	Int id	Model

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	15
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

## 3.2.7 Class-07: DoctorController




## Description:

The DoctorController class is the controller to separates the particular user on role type of doctor and lead to the specific index page about the User model.

No.	Name	Description	Parameter Lists	Return Type
1	findUserModel(Int id)	<ul> <li>id is user id from user table.</li> <li>This method is used to finding the model based on its primary key value of User model.</li> <li>If the model is not found, a 404 HTTP exception will be thrown.</li> <li>This method returns model object by id.</li> </ul>	Int id	Model

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	16
Document Type	Software Design	Release Date	July 30th, 2016	Print Date	July 30 <sup>th</sup> , 2016

## 3.2.8 Class-08: FdaController



## Description:

The FdaController class is the controller to separates the particular user on role type of FDA and lead to the specific index page about the User model.

No.	Name	Description	Parameter Lists	Return Type
1	findUserModel(Int id)	<ul> <li>id is user id from user table.</li> <li>This method is used to finding the model based on its primary key value of User model.</li> <li>If the model is not found, a 404 HTTP exception will be thrown.</li> <li>This method returns model object by id.</li> </ul>	Int id	Model

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	17
Document Type	Software Design	Release Date	July 30th, 2016	Print Date	July 30 <sup>th</sup> , 2016

## **3.2.9 Class-09: Patient**

Patient	
-Spersonal_id : int	
-\$name : string	
-\$surname : string	
-Sprescription : text	
-Screate_by : int	
-\$usemame : string	
Spassword : string	
-Screate_at : datetime	
-\$update_at : datetime	
-\$update_by : string	
+validateIdCard(): Void	

## Description:

This is the Patient model class for table "tbl\_patient". It used to store the patient's information. This model class manages and invokes the patient's information.

#### Attributes:

No.	Name	Description	Type
1	\$personal_id	This attribute is personal ID of patient.	Int
2	\$name	This attribute is name of patients.	String
3	\$surname	This attribute is surname of patients.	String
4	\$prescription	This attribute is prescription.	String
5	\$create_by	This attribute is name of doctor who that create the prescription.	String
6	\$username	This attribute is username of patient.	String
7	\$password	This attribute is password of patient.	String
8	\$create_at	This attribute is date time of prescription that created.	Date Time
9	\$update_by	This attribute is name of doctor who that update the prescription.	String
10	\$update_at	This attribute is date time of prescription that updated.	Date Time

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	18
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

#### Methods:

No.	Name	Description	Parameter Lists	Return Type
1	validateIdCard()	- This method is using for validate the personal id. It invoked in form validator via ajax validation.	-	Void

## 3.2.10 Class-10: PatientController

PatientController	
+getPatientList(): ArrayList	
+viewPatient(Int personal_id) : Mixed	d
+createPatient(Patient patient) : Boo	lean
+updatePatient(Patient patient) : Boo	olean
+deletePatient(Int personal_id) : Boo	olean
+searchPatient(String patients): Arra	ayList
+findModel(Int personal_id) : Model	

## Description:

This class use for controls and manages (CRUD) the information of the patient. This class controls the patient's information by doctor.

No.	Name	Description	Parameter Lists	Return Type
1	getPatientList()	<ul><li>This method is used to providing the list all patients.</li><li>This method is used to display over all list of patients that exist on the database.</li></ul>	-	ArrayList <patient></patient>
2	viewPatient(Int personal_id)	<ul> <li>personal_id is personal id of patient from patient table.</li> <li>This method used for view the single user's information which includes personal id, name, surname, prescription, responsible, username, password and create date to display on the 'view' page by using personal id of patient.</li> </ul>	Int personal_id	Mixed

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	19
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

		<ul> <li>The mechanism of this method involves other method: findModel(int personal_id).</li> <li>This method returns patient's information when personal id that exist on the patient table via user model.</li> <li>This method returns false if the personal id that does not exist on the patient table.</li> </ul>		
3	createPatient(Patien t patient)	- This method is used to creating the new data of patient into database, the data that input to create will be validated by AJAX validation.  - The mechanism of this method	Patient patient	Boolean
		involves other method: validate() and save() from framework.  - This method returns true if doctor input overall information is correct.  -This method returns false if doctor		
		does not input some information or inputting is incorrect.		
4	updatePatient(Patie nt patient)	- This method is used to updating the data of patient from the database by using personal id, the data that input will be validated by AJAX validation.	Patient patient	Boolean
		- The mechanism of this method involves other method: fineModel(int personal_id), validate(), and save() from framework.		
		- This method returns true if doctor input the new overall information is correct.		
		-This method returns false if doctor does not input some information or inputting is incorrect.		

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	20
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

5	deletePatient(Int personal_id)	- personal_id is personal id of patient from patient table.	Int personal_id	Boolean
		- This method is used to deleting the patient from the database by using delete() method from framework via user id.		
		- The mechanism of this method involves other method: findModel(int personal_id) and delete() from framework.		
		- This method returns true when delete personal id that exist on the user table.		
		- This method will return false when delete personal id that does not exist on the user table.		
6	searchPatient(Strin g patients)	- This method is search query applied.	String patients	ArrayList <patient></patient>
		- This method to get ActiveDataProvider instance which will filter models according to data in model fields.		
		- This method returns data items currently available when input attributes that exist in the database.		
6	findModel(Int personal_id)	- personal_id is personal id of patient from patient table.	Int personal_id	Model
		- This method is used to finding the model based on its primary key value.		
		- If the model is not found, a 404 HTTP exception will be thrown.		
		- This method returns model object by personal id.		

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	21
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

## **3.2.11 Class-11: Allergy**

Allergy	
-\$allergy_id : int	
-Spersonal_id : int	
-\$name : string	
-\$surname : string	
-Sreporter : int	
-Screate_at : datetime	
+validateIdCard(): Void	

## Description:

This is the Allergy model class for table "tbl\_allergy". It is using to store the allergy report's information. This model class manages and invokes the allergy report.

#### Attributes:

No.	Name	Description	Type
1	\$allergy_id	This attribute is id of allergy report	Int
2	\$personal_id	This attribute is personal id of the patient.	Int
3	\$name	This attribute is name of the patient.	String
4	\$surname	This attribute is surname of patients.	String
5	\$allergy_drug	This attribute is information of allergy report	String
5	\$reporter	This attribute is name of doctor who are reporter.	String
6	\$create_at	This attribute is date time of report that doctor created.	Date Time

No.	Name	Description	Parameter Lists	Return Type
1	validateIdCard()	- This method is using for validate the personal id. It invoked in form validator via Ajax validation.	-	Void

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	22
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

## 3.2.12 Class-12: AllergyController

ľ	AllergyController
+ç	getAllergyList(): ArrayList
4٧	riewAllergy(Int allergy_id) : Mixed
+ç	getAllergyList2() : ArrayList
+1	riewAllergy2(Int allergy_id) : Mixed
+0	reateAllergy(Allergy allergy) : Boolean
+1	ipdateAllergy(Allergy allergy): Boolean
+0	deleteAllergy(Int allergy_id) : Boolean
+8	searchAllergy(String allergies): ArrayList
	indModel(Int allergy id) : Model

## Description:

This class use for controls and manages (CRUD) the information of the allergy and some patients. This class controls the allergy report's information by doctor.

No.	Name	Description	Parameter Lists	Return Type
1	getAllergyList()	<ul> <li>This method is used to providing the list all allergy reports.</li> <li>This method is used to display over all list of report that exist on the database.</li> <li>This method using for role type doctor.</li> </ul>	-	ArrayList <allergy></allergy>
2	viewAllergy(Int allergy_id)	<ul> <li>- allergy_id is allergy id of allergy report from allergy table.</li> <li>- This method used for view the single allergy report which includes allergy id, personal id, name, surname, allergy drug, reporter and create date to display on the 'view' page by using allergy id of report.</li> <li>- The mechanism of this method involves other method: findModel(int id).</li> <li>- This method returns allergy report when allergy id that exist on the allergy table via user model.</li> </ul>	Int allergy_id	Mixed

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	23
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

3	getAllergyList2()	<ul> <li>This method returns false if the allergy id that does not exist on the allergy table.</li> <li>This method is using for role type doctor.</li> <li>This method is used to providing the list all allergy reports.</li> <li>This method is used to display over all list of report that exist on the database.</li> <li>This method using for role type FDA.</li> </ul>	-	ArrayList <allergy></allergy>
4	viewAllergy2(Int allergy_id)	<ul> <li>allergy_id is allergy id of allergy report from allergy table.</li> <li>This method used for view the single allergy report which includes allergy id, personal id, name, surname, allergy drug, reporter and create date to display on the 'view' page by using allergy id of report.</li> <li>The mechanism of this method involves other method: findModel(int id).</li> <li>This method returns allergy report when allergy id that exist on the allergy table.</li> <li>This method returns false if the allergy id that does not exist on the allergy table.</li> <li>This method is using for role type FDA.</li> </ul>	Int allergy_id	Mixed

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	24
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

5	createAllergy(Aller gy allergy)	<ul> <li>This method is used to creating the new data of allergy report into database, the data that input to create will be validated by AJAX validation.</li> <li>The mechanism of this method involves other method: validate() and save() from framework.</li> <li>This method returns true if doctor input overall information is correct.</li> <li>This method returns false if doctor</li> </ul>	Allergy allergy	Boolean
		does not input some information or inputting is incorrect.		
6	updateAllergy(Alle rgy allergy)	- This method is used to updating the data of allergy report from the database by using allergy id, the data that input will be validated by AJAX validation.	Allergy allergy	Boolean
		- The mechanism of this method involves other method: fineModel(int allergy_id), validate(), and save() from framework.		
		- This method returns true if doctor input the new overall information is correct.		
		-This method returns false if doctor does not input some information or inputting is incorrect.		

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	25
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

7	deleteAllergy (Int allergy_id)	- allergy_id is allergy id of allergy report from allergy table.	Int allergy_id	Boolean
		- This method is used to deleting the allergy report from the database by using delete() method from framework via user id.		
		- The mechanism of this method involves other method: findModel(int allergy_id) and delete() from framework.		
		- This method returns true when delete allergy id that exist on the allergy table.		
		- This method will return false when delete allergy id that does not exist on the allergy table.		
8	searchAllergy(Strin g allergies)	- This method is search query applied.	String allergies	ArrayList
		- This method to get ActiveDataProvider instance which will filter models according to data in model fields.		
		- This method returns data items currently available when input attributes that exist in the database.		
9	findModel(Int allergy_id)	- allergy_id is allergy id of allergy report from allergy table.	Int allergy_id	Model
		- This method is used to finding the model based on its primary key value.		
		- If the model is not found, a 404 HTTP exception will be thrown.		
		- This method returns model object by allergy_id.		

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	26
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

## (Mobile application)

## 3.2.13 Class-13: MainActivity

MainActivity	
toolbar : Toolbar	
onCreate(Bundle savedInstanceState): void onCreateOptionsMenu(enu menu): boolean initInstances(): void	

## Description:

This class is container uses for control and manages the MainFragment class.

#### Attributes:

No.	Name	Description	Type
1	toolbar	Display a toolbar of the main page.	Toolbar

No.	Name	Description	Parameter Lists	Return Type
1	onCreate(bundle savedInstanceState)	This method is use for re-creating activity on the MainFragment called when the activity is starting.	Bundle savedInstance State	void
2	onCreateOptionsMenu()	This method is use to callback method after user click menu on the MainFragment.  - This method returns true when user select button on menu bar.  - This method will return false when user do not select button on menu bar.	Menu menu	boolean
3	initlnstances()	This method is use to handle for click the button on the MainFragment.	-	void

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	27
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

## 3.3.14 Class-14: QrEncodeActivity

QrEncodeActivity
-tvPrescription : TextView
-btnConfirm : Button
-content : String
-tvID : TextView
-tvTime : TextView
-id : String
-prescription : String
-time : String
-REQUEST_CODE : int
+onCreate(Bundle savedInstanceState): void
+initInstances(): void
+intentScanner(): void
+onActivityResult(int requestCode, int resultCode, Intent data): void
+onSaveInstanceState(Bundle state) : void
+onRestoreInstanceState(Bundle state): void
+updateTime(String personalld) : void

## Description:

This class is uses for control the scan QR code function and manage the update time function by pharmacist side.

#### Attributes:

No.	Name	Description	Type
1	tvPrescriptionDetail	Display the text editor of the prescription detail.	TextView
2	tvId	Display the text editor of the ID	TextView
3	tvTime	Display the text editor of the time	TextView
4	btnConfirm	Display the button for confirms.	Button
5	content	Stored the value of content.	String
6	id	Stored the value of id.	String
7	prescription	Stored the value of prescription.	String
8	time	Stored the value of time.	String

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	28
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

No.	Name	Description	Parameter Lists	Retur n Type
1	onCreate(Bundle savedInstanceState)	This method is use for recreating activity called when the activity is starting.	Bundle savedInstanceState	void
2	initInstances()	This method is use to handle for click the button.	-	void
3	intentScanner()	This operation is use for call to use scanner activity.	-	void
4	onActivityResult()	This method is use for receive the requestCode call when activity launched exits.	Int requestCode Int resultCode Intent data	void
5	onSaveInstanceState()	This method is retrieve per- instance state.	Bundle state	void
6	onRestoreInstanceState()	This method is use for receive the SaveInstanceState	Bundle state	void
7	updateTime()	This method is updates the time value to the personal ID column with.	String personalId	void

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	29
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

## 3.3.15 Class-15: LoginActivity

LoginActivity
+onCreate(Bundle savedInstanceState) : void

## Description:

This class is container uses for control and manages LoginFragment class.

No.	Name	Description	Parameter Lists	Return Type
1	onCreate(Bundle savedInstanceState)	This method is use for re-creating activity on the LoginFragment called when the activity is starting.	Bundle savedInstance State	void

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	30
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

## 3.3.16 Class-16: QRCodeEncoder

-WHITE: int
-BLACK: int
-WIDTH: int
-HEIGHT: int
-QUIET\_ZONE\_SIZE: int
-qrlmg: lmageView
+encodeAsBitmap(String str): bitmap
+encode(String contents, BarcodeFormat, int width, int height): BitMatrix
+renderResult(QRCode code, int width, int height, int quietZone): BitMatrix

#### Description:

This class is model use for controls and manage the generate the QR code function of web service.

#### Attributes:

No.	Name	Description	Type
1	WHITE	This attribute use for stored the amount of white box on a QR code.	Int
2	BLACK	This attribute use for stored the amount of black box on a QR code.	Int
3	WIDTH	This attribute use for stored the amount of width on a QR code.	Int
4	HEIGHT	This attribute use for stored the amount of height on a QR code.	Int
5	QUIET_ZONE_SIZE	This attribute use for stored the amount of size on a QR code.	Int
6	qrImg	This attribute use for stored the image of QR code.	ImageView

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	31
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

No.	Name	Description	Parameter Lists	Return Type
1	encodeAsBitmap()	This method use for create the bitmap.	String str	bitmap
2	encode()	This method use for encode a QR code.	String contents, BarcodeFormat format, Int width, Int height	BitMatrrix
3	renderResult()	This method is render the QR code.	QRCode code, Int width, Int height, Int quietZone	BitMatrrix

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	32
<b>Document Type</b>	Software Design	Release Date	July 30th, 2016	Print Date	July 30th, 2016

## 3.3.17 Class-17: ApiService

	ApiService
-BASE_URL : string	
-URL_LOGIN : string	

-URL\_UPDATE\_TIME : string

+post(Context context, String url, RequestParams params, AsyncHttpResponseHandler responseHandler) : void

## **Description:**

This class use for connect with the API service.

#### **Attributes:**

No.	Name	Description	Type
1	BASE_URL	This attribute use for connect the URL of IPv4 address.	String
2	URL_LOGIN	This attributes use for connect login.php of web service.	String
3	URL_UPDATE_TIME	This attributes use for connect update_prescription.php of web service.	String

No.	Name	Description	Parameter Lists	Return Type
1	post()	This method use to request the server.	Context context,  String url, RequestParams params,  AsyncHttpRespons eHandler responseHandler	void

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	33
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

## 3.3.18 Class-18: LoginListener

LoginListener
+onSuccess(String message) : void
+onInternetDown() : void
+onFailed(int statusCode, String responseBody) : void

## Description:

This class is pattern use when the application connects with the login API service.

No.	Name	Description	Parameter Lists	Retur n Type
1	onSuccess()	This method is use for create a message call back to user after the activity successful.	String message	void
2	onInternetDown()	This method is call back to user after cannot the internet.	-	void
3	onFailed()	This method is use for create a message call back to user after the activity fail.	Int statusCode, String responseBody	void

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	34
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

## 3.3.19 Class-19 LoginEventManager

LoginEventManager	
-context : Context	
-manager : PreferencesManager	
+login(final String username, final String password) : void	

## Description:

This class is model use for manage the login function connect with login API service.

## Attributes:

No.	Name	Description	Туре
1	context	This attribute use for stored a data from API.	Context
2	manager	This attribute use for stored a PreferencesManager context.	PreferencesManager

No.	Name	Description	Parameter Lists	Return Type
2	Login()	This method use for login to	String username,	void
		the system.	String password,	

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	35
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

## 3.3.20 Class-20 PreferencesManager

PreferencesManager

-context : Context

-preferences : SharedPreferences -editor : SharedPreferencesEditor

-KEY\_PREFS : string

-KEY\_PERSONAL\_ID : string

-KEY\_FNAME : string -KEY\_LNAME : string

-KEY\_PRESCRIPTION : string

-KEY\_USED : string -KEY\_DATE : string

+login(UserInformationModel userInfo): void

+logout(): void

+getUserInformation(): UserInformationModel

#### Description:

This class uses for stored the information of prescription.

#### Attributes:

No.	Name	Description	Туре
1	context	This attributes use for stored context.	Context
2	preference	This attributes use for stored preference.	SharedPreferences
3	editor	This attributes use for stored editor.	SharedPreferences.Editor
4	KEY_PREFS	This attributes use for stored mPrefs.	String
5	KEY_PERSONAL_ID	This attributes use for stored personal ID.	String
6	KEY_FNAME	This attributes use for stored name.	String
7	EY_LNAME	This attributes use for stored surname.	String
8	KEY_PRESCRIPTION	This attributes use for stored prescription.	String
9	KEY_USED	This attributes use for stored used	String
10	KEY_DATE	This attributes use for stored date.	String

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	36
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

No.	Name	Description	Parameter Lists	Return Type
1	login(UserInformationM odel userInfo)	This method uses for store the information of patient's profile after login successful.	UserInformationM odel userInfo	void
2	logout()	This method uses for clear the value when user logout.	-	void
3	getUserInformation()	This method uses for get the value of UserInformationModel.	-	UserInfor mationM odel

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	37
Document Type	Software Design	Release Date	July 30th, 2016	Print Date	July 30th, 2016

# 3.3.21 Class-21: PrescriptionManager

PrescriptionManager			
-context : Context			
+UpdateTime(final String personalld, final String currentTime) : void			

# Description:

This class is model uses for connect with update\_prescription API service.

# Attributes:

No.	Name	Description	Туре
1	context	This attributes use for stored context.	Context

No.	Name	Description	Parameter Lists	Return Type
1	updateTime()	This method use for update current time to the prescription.	String personalId, String currenTtime	void

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	38
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

# 3.3.22 Class-22: UpdateTimeListener

UpdateTimelistener
ess(String lastTime) : void
etDown()
d(int statusCode, String responseBody)

# Description:

This class is pattern used when connect with the update\_prescription API service.

No.	Name	Description	Parameter Lists	Return Type
1	onInternetDown()	This method is interface for callback invocation on an application action when the internet failed.	-	void
2	onSuccess()	This method is interface for callback invocation on an application action when the operation failed.	String lastTime	void
3	onFailed()	This method is interface for callback invocation on an application action when the operation succeeded	Int statusCode, String responseBody	void

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	39
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### 3.3.23 Class-23: LoginFragment

-context: Contex
-etUsername: EditText
-etPassword: EditText
-btnLogin: Button
-btnScan: TextView
+onCreate(Bundle savedInstanceState): void
+onCreateView(LayoutInflater inflater, ViewGroup container, Bundle savedInstanceState): View
+init(Bundle savedInstanceState)
+initInstances(View rootView, Bundle savedInstanceState): void
+onRestoreInstanceState(Bundle savedInstanceState): void
+onClick(View v): void

### Description:

This class use for controls and manage interface of login page.

#### Attributes:

No.	Name	Description	Туре
1	context	This attributes use for stored context.	Context
2	etUsername	This attributes use for stored username.	EditText
3	etPassword	This attributes use for stored password.	EditText
4	btnLogin	This attributes use for display login button.	Button
5	btnScan	This attributes use for display link scan code button.	TextView

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	40
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

No.	Name	Description	Parameter Lists	Return Type
1	onCreate()	This method called to do initial creation of a LoginFragment.	Bundle savedInstanceState	void
2	onCreateView()	This method called to have the LoginFragment instantiate its user interface view.	LayoutInflater inflater, ViewGroup container, Bundle savedInstanceState	view
3	init()	This method use for set menu item to show in action bar.	-	void
4	initInstances()	This method uses for set the component.	-	void
5	onRestoreInstanceState()	This method uses for restore the value from state.	-	void
6	onClick()	This method is call when a view have been clicked	View v	void

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	41
Document Type	Software Design	Release Date	July 30th, 2016	Print Date	July 30 <sup>th</sup> , 2016

# 3.3.24 Class-24: MainFragment

MainFragment	
context : Context	
manager : PreferencesManager	
userInfo : UserInformationModel	
vName : TextView	
vDate : TextView	
vQrCode : ImageView	
otnQr : Button	
userID : string	
prescription: string	
used : string	
name : string	
name : string	
date: string	
onCreate(Bundle savedInstanceState) : void	
onCreateView(LayoutInflater inflater, ViewGroup container, Bundle savedInstanceStal	te): viev
init(Bundle savedInstanceState) : void	20
initInstances(View rootView, Bundle savedInstanceState) : void	
onRestoreInstanceState(Bundle savedInstanceState) : void	
onClick(View v) : void	
onOptionsItemSelected(MenuItem item) : boolean	

# Description:

This class use for controls and manage interface of main page.

# Attributes:

No.	Name	Description	Туре
1	context	This attribute use for stored a context.	Context
2	manager	This attribute use for stored a manager.	PreferencesManager
3	btnQr	This attribute use for control a button	Button
4	TextView	This attribute use for stored a TextView.	UserInformationModel
5	tvName	This attribute use for display a tvName.	TextView
6	tvDate	This attribute use for display a tvDate.	TextView
7	ivQrCode	This attributes use for display a QrCode.	ImageView
8	userID	This attributes use for stored personal ID.	String
9	prescription	This attributes use for stored prescription.	String
10	used	This attribute use for stored used.	String

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	42
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

11	fname	This attributes use for stored name.	String
12	lname	This attribute use for stored surname.	String

No.	Name	Description	Parameter Lists	Return Type
1	onCreate()	This method use for create a new activity when the application is starting.	Bundle savedInstanceState	void
2	onCreateView()	This method use for define the layout of MainFragment.	LayoutInflater inflater, ViewGroup container, Bundle savedInstanceState	view
3	init()	This method use for set menu item to show in action bar.	-	void
4	initInstances()	This method uses for set the component.	-	void
5	onRestoreInstanceState()	This method uses for restore the value from state.	Bundle savedInstanceState	void
6	onClick()	This method is call when a view have been clicked	View v	void
7	onOptionsItemSelected()	This method uses for control logout system.  - This method returns true when user select logout button.  - This method will return false when user do not select logout button.	MenuItem item	boolean

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	43
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

### 3.3.25 Class-25: Config (web service)



### Description:

This class is web service use for connector between ApiService class and database.


### 3.3.26 Class-26: login (web service)



### Description:

This class is web service use for connector between LoginEventManager class and database.

### **3.3.27** Class-27: update\_prescription (web service)




### Description:

This class is web service use for connector between PrescriptionManager class and database.

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	44
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

# **Chapter Four | Database Design**


# 4.1 ER Diagram



Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	45
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

# **4.2 Table Description**

### 4.2.1 User Table



No.	Field Name	Туре	Length	PK/FK
1	id	int	10	PK
2	name	varchar	45	-
3	surname	varchar	45	-
4	username	varchar	45	-
5	password	varchar	45	-
6	password2	varchar	45	-
7	roleId	int	2	FK


### 4.2.2 Role Table



No.	Field Name	Туре	Length	PK/FK
1	id	tinyint	1	PK
2	name	varchar	30	-

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	46
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016


### 4.2.3 Patient Table



No.	Field Name	Type	Length	PK/FK
1	personal_id	int	17	PK
2	name	varchar	45	-
3	surname	varchar	45	-
4	prescription	longtext	-	-
5	create_by	varchar	45	-
4	username	varchar	45	-
6	password	varchar	45	-
7	create_at	datetime	-	-
8	update_at	datetime	-	-
9	update_by	varchar	45	-
10	userId	int	10	FK

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	47
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

# 4.2.4 Allergy Table



No.	Field Name	Туре	Length	PK/FK
1	allergy_id	int	10	PK
2	personal_id	int	17	-
3	name	varchar	45	-
4	surname	varchar	45	-
5	username	varchar	45	-
6	allergy_drug	longtext	ı	-
7	reporter	varchar	45	-
8	create_at	datetime	1	-
9	update_at	datetime	-	-
10	update_by	varchar	45	-
11	userId	int	10	FK

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	48
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

# **Chapter Five | Sequence Diagram**

SD-01: Doctors can view a doctor home page on the web application.




Figure 1: SD-01 Doctors can view a doctor home page on the web application.

### SD-02: Doctors can create a patient's profile on the web application.




Figure 2: SD-02 Doctors can create a patient's profile on the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	50
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

### SD-03: Doctors can update a patient's profile on the web application.




Figure 3: SD-03 Doctors can update a patient's profile on the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	51
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

### SD-04: Doctors can delete a patient's profile on the web application.




Figure 4: SD-04 Doctors can delete a patient's profile on the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	52
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

### SD-05: Doctors can search a patient's profile on the web application.




Figure 5: SD-05 Doctors can search a patient's profile on the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	53
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

### SD-06: Doctors can view a patient's profile on the web application.




Figure 6: SD-06 Doctors can view a patient's profile on the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	54
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

# SD-07: Doctors can view a list of patient's profile on the web application.




Figure 7: SD-07 Doctors can view a list of patient's profile on the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	55
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

### SD-08: Administrators can view an admin home page on the web application.




Figure 8: SD-08 Administrators can view an admin home page on the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	56
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

### SD-09: Administrators can create a user's profile on the web application.




Figure 9: SD-09 Administrators can create a user's profile on the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	57
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

### SD-10: Administrators can update a user's profile on the web application.




Figure 10: SD-10 Administrators can update a user's profile on the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	58
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

### SD-11: Administrators can delete a user's profile on the web application.




Figure 11: SD-11 Administrators can delete a user's profile on the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	59
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

### SD-12: Administrators can search a user's profile on the web application.




Figure 12: SD-12 Administrators can search a user's profile on the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	60
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

# SD-13: Administrators can view a user's profile on the web application.




Figure 13: SD-13 Administrators can view a user's profile on the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	61
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

# SD-14: Administrators can view a list of user's profile on the web application.




Figure 14: SD-14 Administrators can view a list of user's profile on the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	62
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

### SD-15: Doctors, Administrators, and FDAs can login to the web application.




Figure 15: SD-15 Doctors, Administrators, and FDAs can login to the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	63
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

### SD-16: Doctors, Administrators, and FDAs can logout from the web application.




Figure 16: SD-16 Doctors, Administrators, and FDAs can logout from the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	64
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

### SD-17: Patients can login to the mobile application.




Figure 17: SD-17 Patients can login to the mobile application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	65
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

# SD-18: Patients can logout from the mobile application.




Figure 18: Patients can logout from the mobile application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	66
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

# SD-19: Patients can view the QR code on the mobile application.




Figure 19: SD-19 Patients can view the QR code on the mobile application

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	67
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

# SD-20: Pharmacists can scan the QR code on the mobile application.




Figure 20: SD-20 Pharmacists can scan the QR code on the mobile application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	68
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

### SD-21: Pharmacists can add the time of dispensation to the patient's profile on the mobile application.




Figure 21: SD-21 Pharmacists can add the time of dispensation to the patient's profile on the mobile application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	69
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### SD-22: Doctors can create an allergy report on the web application.




Figure 22: SD-22 Doctors can create an allergy report on the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	70
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

### SD-23: Doctors can update an allergy report on the web application.




Figure 23: SD-23 Doctors can update an allergy report on the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	71
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

### SD-24: Doctors can delete an allergy report on the web application.




Figure 24: SD-24 Doctors can delete an allergy report on the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	72
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

### SD-25: Doctors can search an allergy report on the web application.




Figure 25: of SD-25 Doctors can search an allergy report on the web application.

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	73
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

### SD-26: Doctors can view an allergy report on the web application.




Figure 26: SD-26 Doctors can view an allergy report on the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	74
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### SD-27: Doctors can view a list of allergy report on the web application.




Figure 27: SD-27 Doctors can view a list of allergy report on the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	75
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

### SD-28: FDAs can view a FDA home page on the web application.




Figure 28: SD-28 FDAs can view a FDA home page on the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	76
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

### SD-29: FDAs can view a list of allergy report on the web application.




Figure 29: SD-29 FDAs can view a list of allergy report on the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	77
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

### SD-30: FDAs can view a list of allergy report on the web application.




Figure 30: SD-30 FDAs can view a list of allergy report on the web application.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	78
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

# **Chapter Six | User Interface Design**

**UI-01: Main page** 

Smart Prescription Login

# Welcome!

Smart Prescription Application


© Smart Prescription Application 2016

This page is main page use for administration, doctor, and FDA that consist of

- (1) "Login" button use for login to the web application.
- (2) "Home" button use for moving to the Home page.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	79
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### **UI-02:** Administration home page



© Smart Prescription Application 2016

This page is administration home page for administration that consist of

- (1) "Logout" button use for logout from the web application.
- (2) "Home" button use for moving to the Home page.
- (3) Information of user that working on the web application.
- (4) Option menu for create users and manage users.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	80
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### **UI-03:** Create user page

Smart Prescription		Home Logout (admin)
Home / Manage Users / Create User		
Create User		
Name		
Surname		
Username		
Password		
Password2		
Position		
	Create	


© Smart Prescription Application 2016

This page is create user page for administration that consist of

- (1) "Logout" button use for logout from the web application.
- (2) "Home" button use for moving to the Home page.
- (3) Breadcrumb link use for select to the previous page.
- (4) Users create form for input the information of user's profile.
- (5) "Create" button for select to create user's profile.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	81
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### **UI-04: Update user page**




This page is update user page for administration that consist of

- (1) "Logout" button use for logout from the web application.
- (2) "Home" button use for moving to the Home page.
- (3) Breadcrumb link use for select to the previous page.
- (4) Users update form for input the information of user's profile.
- (5) "Update" button for select to update user's profile.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	82
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### **UI-05:** Manage users page



This page is manage user page for administration that consist of

- (1) "Logout" button use for logout from the web application.
- (2) "Home" button use for moving to the Home page.
- (3) Breadcrumb link use for select to the previous page.
- (4) The number and total amount of user's profiles.
- (5) List of user's profile include six columns (no., position, name, surname, username, and options.).
- (6) "Create User" button for select to create user's profile.
- (7) "Options column of table" uses for view, edit, and delete the user's profile.
- (8) Search text field use for input a name, surname or username of user's profile.
- (9) "Search" button for select to search user's profiles.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	83
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### **UI-06:** User's profile page

# Home / Manage Users / Dr.Jazz Dr. Jazz ID 32 Name Dr.Jazz Surname Hamonica Username jazz Password jazz Position Doctor




This page is user's profile page for administration that consist of

- (1) "Logout" button use for logout from the web application.
- (2) "Home" button use for moving to the Home page.
- (3) Breadcrumb link use for select to the previous page.
- (4) Information of user's profile.
- (5) "Update" button use for select to update the user's profile.
- (6) "Delete" button use for select to delete the user's profile.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	84
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### **UI-07: Doctor home page**




This page is doctor home page for doctor that consist of

- (1) "Logout" button use for logout from the web application.
- (2) "Home" button use for moving to the Home page.
- (3) Information of doctor that working on the web application.
- (4) Option menu for create patient and manage patient.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	85
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### **UI-08:** Create patient page




This page is create patient page for doctor that consist of

- (1) "Logout" button use for logout from the web application.
- (2) "Home" button use for moving to the Home page.
- (3) Breadcrumb link use for select to the previous page.
- (4) Patients create form for input the information of patient's profile.
- (5) "Create" button for select to create a patient's profile.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	86
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### **UI-09:** Update patient page




This page is update patient page for doctor that consist of

- (1) "Logout" button use for logout from the web application.
- (2) "Home" button use for moving to the Home page.
- (3) Breadcrumb link use for select to the previous page.
- (4) Patient update form for input the information of patient's profile.
- (5) "Update" button for select to update patient's profile.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	87
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### **UI-10:** Manage patient page




This page is manage patient page for doctor that consist of

- (1) "Logout" button use for logout from the web application.
- (2) "Home" button use for moving to the Home page.
- (3) Breadcrumb link use for select to the previous page.
- (4) The number and total amount of patient's profiles.
- (5) List of patient's profile include seven columns (no., personal ID, name, surname, responsible, create date, and options.).
- (6) "Create patient" button for select to create patient's profile.
- (7) "Options column of table" uses for view, edit, and delete the patient's profile.
- (8) Search text field use for input a name, surname or personal ID of patient's profile.
- (9) "Search" button for select to search patient's profiles.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	88
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### **UI-11: Patient's profile page**



This page is patient's profile page for doctor that consist of

- (1) "Logout" button use for logout from the web application.
- (2) "Home" button use for moving to the Home page.
- (3) Breadcrumb link use for select to the previous page.
- (4) Information of patient's profile.
- (5) "Update" button use for select to update the patient's profile.
- (6) "Delete" button use for select to delete the patient's profile.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	89
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### **UI-12:** Create allergy report page


Smart Prescription		Home	Logout (doctor)
Home / Allergies / Create Allergy			
Create Allergy			
Personal ID			
Name			
Surname			
Allergy Drug			
	Create		

This page is create allergy report page for doctor that consist of

- (1) "Logout" button use for logout from the web application.
- (2) "Home" button use for moving to the home page.
- (3) Breadcrumb link use for select to the previous page.
- (4) Allergy reports create form for input the information of allergy report.
- (5) "Create" button for select to create allergy report.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	90
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### **UI-13:** Update allergy report page




This page is update allergy report page for doctor that consist of

- (1) "Logout" button use for logout from the web application.
- (2) "Home" button use for moving to the Home page.
- (3) Breadcrumb link use for select to the previous page.
- (4) Allergy report updates form for input the information of allergy report.
- (5) "Update" button for select to update allergy report.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	91
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### **UI-14:** Manage allergy reports page




This page is manage allergy report page for doctor that consist of

- (1) "Logout" button use for logout from the web application.
- (2) "Home" button use for moving to the home page.
- (3) Breadcrumb link use for select to the previous page.
- (4) The number and total amount of allergy report.
- (5) List of allergy report include six columns (no., personal ID, name, surname, reporter, create date, and options.).
- (6) "Create report" button for select to create allergy report.
- (7) "Options column of table" uses for view, edit, and delete the allergy report.
- (8) Search text field use for input a name, surname or personal ID of allergy report.
- (9) "Search" button user for select to search allergy report.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	92
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016

### **UI-15:** Allergy report page




This page is allergy report page for doctor that consist of

- (1) "Logout" button use for logout from the web application.
- (2) "Home" button use for moving to the home page.
- (3) Breadcrumb link use for select to the previous page.
- (4) Information of allergy report.
- (5) "Update" button use for select to update the allergy report.
- (6) "Delete" button use for select to delete the allergy report.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	93
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### **UI-16: FDA home page**




This page is FDA home page for FDA that consist of

- (1) "Logout" button use for logout from the web application.
- (2) "Home" button use for moving to the Home page.
- (3) Information of FDA that working on the web application.
- (4) Option menu for view allergy report.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	94
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### **UI-17: List allergy report page**




This page is list allergy report page for FDA that consist of

- (1) "Logout" button use for logout from the web application.
- (2) "Home" button use for moving to the home page.
- (3) Breadcrumb link use for select to the previous page.
- (4) The number and total amount of allergy report.
- (5) List of allergy report include six columns (no., personal ID, name, surname, reporter, create date, and options.).
- (6) "Options column of table" uses for view the allergy report.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	95
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### **UI-18: FDA Allergy report page**




This page is FDA allergy report page for FDA that consist of

- (1) "Logout" button use for logout from the web application.
- (2) "Home" button use for moving to the home page.
- (3) Breadcrumb link use for select to the previous page.
- (4) Information of allergy report.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	96
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### **UI-19: Login page**




This page is Login page for doctor, administration, and FDA that consist of

- (1) "Login" button use for moving to the login page.
- (2) "Home" button use for moving to the home page.
- (3) Login form for input a username and password.
- (4) "Login" button for click to login.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	97
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### (Mobile Application)

### **UI-20:** Login page (mobile)



This page is Login page (mobile) for patient that consist of

- (1) Login form for input a username and password.
- (2) "Login" button for click to login.
- (3) "Scan QR code" button for click to scan QR code.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	98
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

# **UI-21: Main page (mobile)**



This page is main page (mobile) for patient that consist of

- "Logout" button for logout from the application.
   Name of patient's profile.
- (3) "SHOW QR CODE" button for click to display QR code.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	99
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### **UI-22: QR code page (mobile)**



This page is QR code page for patient that consist of

- (1) "Logout" button for logout from the application.
- (2) Name of patient's profile.(3) "SHOW QR CODE" button for click to display QR code.
- (4) QR code of the patient's profile.

Document Name	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	100
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30th, 2016


# **UI-23: Scan QR code page (mobile)**



This page is scan QR code page for pharmacist that consist of (1) Scan QR code function.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	101
<b>Document Type</b>	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016

### **UI-24: Pharmacist page (mobile)**



- This page is pharmacist page for pharmacist that consist of
  (1) Information of patient's profile.(2) "CONFIRM TO USE" button for click to record the current time.

<b>Document Name</b>	SMART-PRESCRIPTION- APPLICATION-[Software- Design _V2.0]	Owner	Natthakn Kaeokanpai, Phithiwat Sitthitun	Page	102
Document Type	Software Design	Release Date	July 30 <sup>th</sup> , 2016	Print Date	July 30 <sup>th</sup> , 2016